CapsoVision   
*empowering through innovation*

Capso  CAM  
PLUS<sup>®</sup>

## INSTRUCTIONS FOR USE

# CapsoCamPlus<sup>®</sup>

English Deutsch Español Français Italiano Nederlands

CE 0197

## Indications for Use

The CapsoCam Plus® capsule system is intended for visualization of the small bowel mucosa in patients ages 2 and above. It may be used as a tool in the detection of abnormalities of the small bowel.

## Indications for Pediatric Use

Based upon review of the published literature and available guidelines published to date, indications for Capsule Endoscopy (CE) in pediatric patients are as follows:

- Diagnosis of Inflammatory Bowel Disease (IBD), including the evaluation of suspected Crohn's Disease (CD) and the ongoing diagnostic surveillance of pediatric patients with established CD.
  - A normal CE study has a high specificity for excluding small bowel CD.
  - CE should be recommended in suspected cases of CD when conventional endoscopy and imaging tools are not feasible or have been non-diagnostic and no bowel obstructive symptoms are present.
  - The possibility of capsule retention should be carefully considered in patients at risk for intestinal stenosis, mainly those with a previous abdominal surgery and those with suspected or established Crohn's disease.
- Diagnosis of Occult Gastrointestinal Bleeding (OGIB) and/or Chronic Anemia—as an alternative diagnostic methodology and as a tool for management of OGIB in pediatric patients.
  - Diagnosis of Polyposis Syndromes—providing safe, accurate detection of small-bowel polyps in pediatric patients. Specifically, CE may be a useful tool in the screening and surveillance of patients with Peutz-Jeghers syndrome (PJS).
  - CE may be useful as an alternative diagnostic modality for other diseases affecting the pediatric population such as: celiac disease and malabsorption syndromes, lymphangiectasia, graft-versus-host disease, chronic abdominal pain, and failure to thrive.

## Device Description

CapsoCam Plus is a single-use, ingestible capsule that acquires and stores endoscopy examination images in on-board memory while moving through the gastrointestinal tract, propelled by natural peristalsis. The patient retrieves the capsule using the provided retrieval kit and returns it to the physician who downloads and reviews the images on a computer. The capsule is typically excreted within 3 to 30 hours after being swallowed.

## Ratings

IP68

Type BF applied part

3.1V DC battery-powered

10mA

Class IIa

## Device Common Name

System, Imaging, Gastrointestinal, Wireless, Capsule

## Operating Range

Device operates in vivo. Ex vivo it operates within the following environmental ranges—

Temperature: 5°–40° C

Humidity: 5% to 95%

Atmospheric Pressure: 59kPa to 105kPa

## Operating Environment

The intended use environments for CapsoCam Plus are home healthcare environment, hospital, transport and ambulatory environments.

## Contraindications

The CapsoCam Plus capsule is contraindicated in patients:

- Who have known or suspected gastrointestinal obstructions, strictures or fistula
- Who are pregnant
- Who have gastroparesis
- Who have a swallowing disorder

## Patient Condition

The physician should consider performing a small bowel series before utilizing this device in patients who are suspected to have strictures or fistulas. See Risks on Page 4 for additional risks related to specific patient conditions.

## Adverse Events

Potential adverse events associated with the use of this device may include:

- Obstruction, perforation, and mucosal injury or bleeding.
- Delayed or no excretion of the capsule, and aspiration. In some instances, intervention is required to remove the capsule.

## Cautions

- Ensure that only personnel, who are trained on all of the CapsoCam Plus Capsule Endoscope System operating procedures use the system.
- In a small number of cases, the CapsoCam Plus capsules may not image the entire small bowel due to variations in patient GI motility or anatomy.
- Final diagnosis based on the video images should be made only by physicians who are trained in the interpretation of capsule endoscopy images.

## Benefits

- The CapsoCam Plus Capsule Endoscopy System operates without the use of radiofrequency (RF) transmission. All of the diagnostic data is stored on board the capsule, which allows patients to undergo the examination unhindered by external receiver equipment.
- Due to the lack of RF transmission, patients with an implantable device can safely ingest the capsule.
- It allows your doctor to examine the small bowel which cannot be reached by traditional upper endoscopy or colonoscopy.
- It's a simple, safe and less invasive alternative to traditional imaging procedures. The procedure does not require sedation or radiation. Patients may continue with their normal daily activity during the procedure.

## Risks

- All medical procedures carry some risks. Information in this manual should not be used as a substitute for talking with your doctor about diagnosis and treatment.
- CapsoCam Plus capsule endoscopy is not for everyone. CapsoCam Plus is contraindicated in patients with known or suspected gastrointestinal obstruction, strictures or fistulas, patients who are pregnant, patients with gastroparesis and patients with swallowing disorders.

- Capsule retention has been reported in less than 2% of all capsule endoscopy procedures. Capsule retention is defined as having a capsule remain in the digestive tract for more than 72 hours.
- Causes of capsule retention cited in the literature include: NSAID strictures, Crohn's disease, small bowel tumors, intestinal adhesions, ulcerations, and radiation enteritis.
- Although each risk has been addressed and mitigated as far as possible, there is a possibility that the CapsoCam Plus can fail due to components, hardware, software and packaging issues which may cause delay of diagnosis, image quality degradation, incomplete exam, or data corruption.
- The risks of capsule leakage, allergic reaction, biological contamination, and infection have been addressed and mitigated as far as possible. There is a possibility that these can occur if care is not taken in handling the product according to these instructions for use and hospital procedures, but the probability of occurrence is very low.
- Summaries in published literature identify the risk of capsule retention to be approximately 1.5% for obscure bleeding, 1.4% for suspected Crohn's disease, 5% for known Crohn's, and 2.1% for neoplastic lesions. [1], [2].  
[1] Cave et al. *Endoscopy*. 2005; 37: 1065-1067.  
[2] Zhuan et al. *GI Endoscopy*. 2010; 71: 280-286.

## Special Considerations for Pediatrics

- Consideration should be given to the risk versus benefit of performing a small bowel series in pediatric patients to evaluate the presence or absence of strictures or fistulas in the small bowel of pediatric patients since the risk of the potentially detrimental effect of radiation may be greater to the pediatric patient and the diagnostic accuracy of this procedure may not effectively exclude the potential of capsule retention. Magnetic Resonance Endoscopy (MRE) or the administration of a Patency Capsule should be considered as possible diagnostic modalities to evaluate the potential presence of strictures or fistulas and the potential risk of capsule retention.
- Consideration should be given to the risk versus benefit of utilizing a commercially available patency capsule for evaluation of small bowel patency for CapsoCam Plus capsule (30.5mm x 11.3mm) in pediatric patients. While patency capsules are useful tools for evaluating small bowel patency, they do not always fully characterize the presence of possible strictures in the small bowel. Patency capsules may also result in additional risk of infection and perforation as they may be retained in the small bowel. The manufacturer's recommendations regarding pediatric use of commercially available patency capsules should be consulted.
- Consideration should be given to the risk versus benefit of using CE in pediatric patients as the initial diagnostic

modality for evaluation of suspected OGIB and/or chronic anemia because there are no clinical evaluations specifically addressing the indication or performance of CE in children with OGIB. Clinical guidelines published to-date suggest that gastroscopy or colonoscopy be utilized as initial diagnostic modalities in pediatric patients suspected or documented to have OGIB and/or chronic anemia. CE may be effective as a diagnostic modality for use when other methods have failed to document causation of OGIB and/or chronic anemia.

- Consideration should be given to the risk versus benefit of using CE in pediatric patients in various states of malnutrition as there are some limited data from a single source suggesting that pediatric patients with lower BMI may be at increased risk of capsule retention.
- Complications of CE in the pediatric patient, while seldom seen based upon review of published guidelines, are similar to those in adult patients, with capsule retention being the most significant complication. The possibility of capsule retention should be considered when prescribing CE in the pediatric population.
- Additional complications that have been observed in pediatric uses are: nausea, perceived inability to swallow the capsule, delayed transit time of unknown etiology, and temporary gastric retention.

## Bowel Preparation

Patient should follow physician's instructions for bowel preparation.

## Preparation for Use

- Inspect the foil seal. Do not use if any breakage in the foil seal or in the foil itself is observed.
- Prepare a glass or bottle of water. The recommended amount of water is 0.5 liter (~16 fluid ounces).

## Directions for Use

1. Open the capsule package by peeling back the foil cover.
2. Using gloves, remove the plastic lid covering the capsule.
3. Grasp the capsule carefully and then pull it out and away from the package. **Take care not to drop the capsule on the floor.** Dropping the capsule may damage it, which could result in an incomplete examination.
4. Within approximately ten seconds of being removed from the magnet, the capsule LEDs will light up in a consistent rhythm of approximately 5 flashes per second, indicating that the capsule is activated. Because CapsoCam Plus uses auto-illumination to optimize lighting conditions for in vivo images, it is possible that not all capsule LEDs will appear illuminated. If unsure that the capsule is flashing, enclose the capsule in your hand to provide a dark environment to confirm that the LEDs are illuminated.

Once you have confirmed that the capsule is flashing, the capsule is ready to be swallowed.

If the capsule has not begun flashing within 20 seconds after being removed from the package, place the battery-end of the capsule against the magnet in the package and repeat step 4 again. If the capsule still fails to begin flashing, the capsule is malfunctioning and should be replaced.

5. The capsule must be swallowed within 10 minutes of being removed from the package. The capsule should be placed inside the patient's mouth and swallowed with the entire glass of water. Patients should avoid biting the capsule. If the patient is unable to swallow the capsule within 10 minutes, return the capsule to the package with the battery-end of the capsule touching the magnet inside the package, which returns the capsule to the "OFF" state.
6. The patient may engage in normal daily activities while the capsule moves within the digestive tract. The patient should consult the physician if he or she has questions about engaging in a particular activity.
7. Notify the patient that the capsule is typically excreted between 3 and 30 hours after being swallowed. Patient must retrieve the capsule by following the Capsule Retrieval Kit Instructions for Use.

## Warnings

Laser Radiation— Class 1M Laser Do not view directly with optical instruments

- The capsule should not be swallowed by patients for whom a concern of aspiration exists (e.g. due to cognitive or neurological deficits or a history of aspiration). In these patients, it is recommended that a capsule endoscopic delivery system be used to place the capsule directly in the duodenum. Placement of the capsule in the duodenum will decrease the risk of aspiration of the device [by vomiting] and gastric retention.
- A patient with known or suspected delayed gastric emptying (whether disease-related or drug-induced) could be at increased risk for an incomplete examination.
- Instruct the patient to contact the physician immediately if, after ingesting the CapsoCam Plus capsule, there is any abdominal pain, nausea, or vomiting.
- Only one CapsoCam Plus capsule should be ingested at a time and only after confirmation that no other CapsoCam Plus capsules remain in the patient's body.
- If, contrary to instructions, a patient ingests more than one CapsoCam Plus capsule, instruct the patient to immediately contact the physician.
- In patients with unsuspected strictures of the GI tract, any CapsoCam Plus capsule can potentially cause intestinal obstruction resulting in the need for hospitalization and surgery.
- The safety of this device in pregnant women has not been established.
- The safety of this device in patients with significant gastrointestinal diverticular disease is unknown.
- Final diagnosis should be made only by physicians who are trained in the interpretation of capsule endoscopy images.
- The CapsoCam Plus Capsule Endoscopy System needs to be installed and put into service according to the Electromagnetic Compatibility (EMC) information provided in the accompanying documents.
- Patients should avoid biting the CapsoCam Plus capsule prior to swallowing.
- If the capsule is damaged in any way, including by forceful biting, the patient should not swallow the capsule.
- If the capsule is not excreted after 72 hours, patient should contact the physician.

- Store CapsoCam Plus capsules and their packages in a safe place, out of the reach of children and infants.
- If a person has accidentally swallowed a CapsoCam Plus capsule, or the embedded package magnets, seek immediate medical attention.
- Keep package, including embedded magnet, at least 5cm (1.9685 inches) away from pacemakers and other active-implant medical devices.
- Do not modify the capsule without authorization from the manufacturer.
- Patients should not board an aircraft until the capsule has been excreted.
- CapsoCam Plus needs special precautions regarding EMC and needs to be utilized according to the EMC information provided in this manual.
- Portable and mobile RF communications equipment can affect CapsoCam Plus.
- The use of accessories, transducers and cables other than those specified by CapsoVision, Inc. may result in increased EMISSIONS or decreased IMMUNITY of the CapsoCam Plus.
- CapsoCam Plus should not be used adjacent to or stacked with other equipment and if adjacent or stacked use is necessary, CapsoCam Plus should be observed to verify normal operation in the configuration in which it will be used.
- A negative or normal result obtained by the CapsoCam Plus capsule does not exclude the presence of pathology and if symptoms persist, further evaluation should be performed.
- If intestinal fistulas, strictures, or stenoses are suspected, or the patient has had prior abdominal or pelvic surgery, the physician should consider performing a procedure to ascertain patency for an object the size of the CapsoCam Plus capsule.
- After ingesting the CapsoCam Plus capsule and until it is excreted, the patient should not be near any source of powerful electromagnetic fields.
- Patients should not have an MRI exam performed until the capsule has been excreted. Possible patient injury and medical complications could occur.
- The CapsoCam Plus capsule should not be administered to patients below two (2) years of age.
- Thorough understanding of the technical principles, clinical applications and risks associated with the CapsoCam Plus Capsule Endoscopy System is necessary before using this product. Read the entire manual before using the system for the first time.
- To prevent the patient from being exposed to unforeseen risks during passage of the CapsoCam Plus capsule, make sure the patient thoroughly understands the procedure,

and provide the patient with a copy of the Patient Instructions.

- Do not use CapsoCam Plus after its expiration date.
- When swallowing the capsule, there is a possibility of choking on the capsule. If the patient exhibits any symptoms and/or clinical signs of choking (labored breathing, wheezing, involuntary coughing, etc.), the capsule shall be ingested under the supervision of authorized medical personnel.
- Used capsules are biologically contaminated and if reused can cause infection from bacteria, viruses, or other agents.

## Storage

Store the CapsoCam Plus capsule under normal indoor environmental conditions. Do not remove the CapsoCam Plus capsule from the packaging until just prior to use.

## Cleaning

Retrieved capsules should be handled with gloves and cleaned and disinfected as follows:

1. Clean the capsule using ENZOL® Enzymatic Cleaner or equivalent, according to the manufacturer's instructions. Use a soft bristle brush (paying special attention to the dimple at the bottom of the capsule) as necessary to thoroughly remove all debris.
2. Disinfect the capsule using Revital-Ox™ Resert® High Level Disinfectant or equivalent and rinse thoroughly, according to the manufacturer's instructions.

3. Dry capsule completely.

## Data Download at Clinic

Prior to downloading data from the used capsule, it must be cleaned, disinfected, and dried. See the CapsoView® Instructions for Use.

## Disposal

Dispose of capsule and packaging as per local ordinances.

## Warranty

CapsoVision warrants that this capsule is free from defects in both materials and workmanship. Suitability for use of the capsule functionality for any procedure shall be determined by the user. CapsoVision shall not be liable for incidental or consequential damages of any kind. The above warranties are in lieu of all other warranties, either express or implied, including any warranty of any merchantability or fitness for use.

## Essential Performance

Data integrity during image capture and download, laser Class 1M, and single-use are considered to be the essential performance of the CapsoCam Plus capsule.

## Federal Communication Commission (FCC) Compliance

The CAPSOVISION system (CapsoCam Plus) complies with Part 15 of the United States FCC rules and with international standards for electromagnetic compatibility regarding its use.

(Table 1)

<b>Guidance and manufacturer's declaration – electromagnetic emissions</b>		
The CapsoCam Plus is intended for use in the electromagnetic environment specified below. The customer or the user of the CapsoCam Plus should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The CapsoCam Plus uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The CapsoCam Plus is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.

(Table 2)

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The CapsoCam Plus is intended for use in the electromagnetic environment specified below. The customer or the user of the CapsoCam Plus should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
(50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical domestic, commercial or hospital environment.

(Table 3)

<b>Guidance and manufacturer's declaration – electromagnetic immunity</b>			
The CapsoCam Plus is intended for use in the electromagnetic environment specified below. The customer or the user of the CapsoCam Plus should assure that it is used in such an environment.			
<b>Immunity test</b>	<b>IEC 60601 Test Level</b>	<b>Compliance level</b>	<b>Electromagnetic environment – guidance</b>
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not Applicable	Portable and mobile RF communications equipment should be used no closer to any part of the EQUIPMENT, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	20 V/m	<p><b>Recommended separation distance</b></p> <p><math>d = 0.18\sqrt{P}</math> 80 MHz to 800 MHz  <math>d = 0.35\sqrt{P}</math> 800 MHz to 2.5 GHz</p> <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey <sup>a</sup>, should be less than the compliance level in each frequency range.<sup>b</sup></p>

(Table 3 continued)

			<p>Interference may occur in the vicinity of equipment marked with the following symbol:</p> 
<p>NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			
<p><sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the CapsoCam Plus is used exceeds the applicable RF compliance level above, the CapsoCam Plus should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the CapsoCam Plus.</p> <p><sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.</p>			

(Table 4)

**Recommended separation distances between  
portable and mobile RF communications equipment and the CapsoCam Plus**

The CapsoCam Plus is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the CapsoCam Plus can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the CapsoCam Plus as recommended below, according to the maximum output power of the communications equipment.

Rated Maximum output power of transmitter		Separation distance according to frequency of transmitter	
W	m		
	80 MHz to 800 MHz $d = 0.18\sqrt{P}$	800 MHz to 2.5 GHz $d = 0.35\sqrt{P}$	
0.01	0.02	0.04	
0.1	0.06	0.11	
1	0.18	0.35	
10	0.57	1.1	
100	1.8	3.5	

For transmitters rated at a maximum output power not listed above, the recommended separation distance ( $d$ ) in meters ( $m$ ) can be estimated using the equation applicable to the frequency of the transmitter, where ( $P$ ) is the maximum output power rating of the transmitter in watts ( $W$ ) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Symbol	Meaning	Symbol	Meaning
	Manufacturer's catalogue number		Do not reuse this product
	Caution: Consult accompanying documents		Contact/European representative of manufacturer
	Product meets European standards for safety and quality		Use by date (with accompanying date)
	Instructions are included and must be read		Item serial number (with accompanying number)
	Special disposal for electronic waste required		Manufacturer
	Prescription use only For sale by or on the order of a physician		Laser Radiation—Class 1M Laser Product Do not view directly with optical instruments
	Item Lot Number (with accompanying number)		Type BF applied part
	Temperature Limitation		MR unsafe
	Date of manufacture		Quantity
	Atmospheric pressure limitation		Do not use if package is damaged
	Humidity limitation		

#### Manufacturer

CapsoVision, Inc.  
18805 Cox Ave, Suite 250  
Saratoga, CA 95070 USA  
Tel: +1 408 624 1488  
Email:  
customerservice@capsovision.com

#### European Authorized Representative

Dr. Hans-Joachim Lau  
Airport Center (Building C)  
Flughafenstraße 52a  
22335 Hamburg, Germany  
Fax: +49 40 53299 100

#### Conformity Assessment Body

TÜV Rheinland LGA Products GmbH  
90431 Nuremberg, Germany

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